

# ASTOCOVER






## Passive warming insulation

The ASTOCOVER can be used in pre-operative, peri-operative and post-operative situations, in intensive care or on the ward, as prophylaxis against the patient cooling down and as a treatment for hypothermia.

With ASTOCOVER it is prevented that warmth is radiated from the patient into the area and it comes thus to no unnecessary calorific loss of the patient. Hence the patient is protected against cooling down caused by operation theatre ventilation and draught.

An ideal solution to every application can be found in the above-mentioned scenarios by selecting applied parts which can be combined.

In combination with our warming systems such as ASTOPAD SYSTEM ( Patient Warmer ), ASTOFLO and ASTOTHERM PLUS (Blood and Infusion Warmers) a complete conception is realized to prevent the patient from the negative side effects of an operation such as hypothermia and decubitus.

	<b>ASTOCOVER PAS100 Upper body cover</b> to passive warming insulation, 100 X 100 cm <u>Order no.:</u> PAS 100	
	<b>Cover can be washed and disinfected</b> for ASTOCOVER PAS100, 100 X 100 cm <u>Order no.:</u> COV 40100	
	<b>ASTOCOVER PAS150 Body cover</b> to passive warming insulation, 50 X 150 cm <u>Order no.:</u> PAS 150	
	<b>Cover can be washed and disinfected</b> for ASTOCOVER PAS150, 50 X 150 cm <u>Order no.:</u> COV 40150	
	<b>ASTOCOVER PAS155 Arm - Chest cover</b> to passive warming insulation, 50 X 150 cm <u>Order no.:</u> PAS 155	
	<b>Cover can be washed and disinfected</b> for ASTOCOVER PAS150, 50 X 150 cm <u>Order no.:</u> COV 40155	

**Material:** The cover of ASTOCOVER consists of a PU coated PES fabric.

Characteristics:

- breathe-active
- virus and bacteria impermeable
- blood, secretion and urine impermeable
- PVC free
- flexible and no formation of wrinkles
- Ökotex 100, II Certificate

**Cleaning and Disinfection:**

For cleaning and disinfection usual, alcohol-based cleaning agents and disinfectants may be used. The washable and disinfected covers for ASTOCOVER (COV 40100/40150/40155) additionally may be washed at 95°C.



**STIHLER ELECTRONIC**

Stihler Electronic GmbH  
Julius-Hölder-Straße 36 • D-70597 Stuttgart • Germany  
www.stihlerelectronic.de • E-Mail: info@stihlerelectronic.de  
national: Telefon 07 11- 72 06 70 • Telefax 07 11- 720 67 57  
international: Telefon ++ 49- 711- 72 06 70 • Telefax ++ 49- 711- 720 67 57

*Technical changes reserved*