

ASTOPAD[®] System

Heating Elements

ASTOPAD COV070 - Heated blanket

Particularly for application with babies and infants.

Operational areas:

- Nursery
- Delivery room
- Operating theatre

Dimension: 680 X 480 mm

Order No.: COV070

Cover for heated blanket COV070 can be washed (95°C) and disinfected

Order No.: COV40070

Cotton cover for heated blanket COV070 can be washed (95°C)

Order No.: COV20070



ASTOPAD COV105 - Heated blanket

Particularly for application at urological interventions and caesareans.

Dimension: 1050 X 500 mm

Order No.: COV105

Cover for heated blanket COV105 can be washed (95°C) and disinfected

Order No.: COV40105



ASTOPAD COV150 - Heated blanket

Universal use:

- As an underblanket on the operating table
- As an upper body blanket
- As a lower body blanket...

Dimension: 1500 X 500 mm

Order No.: COV150

Cover for heated blanket COV150 can be washed (95°C) and disinfected

Order No.: COV40150

Disposable cover for COV150, VE=50 Stck.

Order No.: COV30150



ASTOPAD COV155 - Arm-Chest-Heated blanket

Heating element to cover the upper body

Dimension: 1500 X 500 mm

Order No.: COV155

Cover for heated blanket COV155 can be washed (95°C) and disinfected

Order No.: COV40155



ASTOPAD COV180 - Heated blanket

Full-body blanket for universal use.

- As a heated blanket under the patient
- and on the top as a top blanket

Dimension: 1800 X 800 mm

Order no.: COV180

Cover for heated blanket COV180 can be washed (95°C) and disinfected

Order No.: COV40180

Disposable cover for COV180, VE=25 Stck.

Order No.: COV30180



ASTOPAD OPT130 - Heatable operating table pad

Is the thoroughly tried and tested heatable operating table pad.

Dimension: 1000 X 500 mm

Order No.: OPT130



STIHLER ELECTRONIC

Stihler Electronic GmbH

Julius-Hölder-Straße 36 • D-70597 Stuttgart • Germany

www.stihlerelectronic.de • E-Mail: info@stihlerelectronic.de

national: Telefon 07 11- 72 06 70 • Telefax 07 11- 720 67 57

international: Telefon ++ 49- 711- 72 06 70 • Telefax ++ 49- 711- 720 67 57

ASTOPAD[®] System, ASTOCOVER and ASTOGEL

Application Recommendations

(COV150, COV180 and PAS155)



ASTOGEL OPG130
Gel pad for prophylaxis of pressure sores
ASTOPAD COV150
Heated blanket as under blanket
ASTOPAD COV150
Heated blanket as arm-chest cover



ASTOGEL OPG130
Gel pad for prophylaxis of pressure sores
ASTOPAD COV150
Heated blanket as under blanket
ASTOPAD COV150
Heated blanket as leg bag



ASTOGEL OPG130
Gel pad for prophylaxis of pressure sores
ASTOPAD COV150
Heated blanket as under blanket
ASTOPAD COV150
Heated blanket as leg bag



ASTOPAD COV150
Heated blanket as leg bag
ASTOCOVER PAS155
Passive warming insulation
as arm-chest cover



ASTOGEL OPG130
Gel pad for prophylaxis of pressure sores
ASTOPAD COV150
Heated blanket as under blanket
ASTOPAD COV150
Heated blanket as leg bag
ASTOCOVER PAS155
Passive warming insulation
as arm-chest cover



ASTOGEL OPG130
Gel pad for prophylaxis of pressure sores
ASTOPAD COV150
Heated blanket as under blanket
ASTOPAD COV150
Heated blanket as leg bag
ASTOCOVER PAS155
Passive warming insulation
as arm-chest cover



ASTOGEL OPG130
Gel pad for prophylaxis of pressure sores
ASTOPAD COV150
Heated blanket as under blanket
ASTOPAD COV180
Heated blanket as full body cover



ASTOGEL OPG130
Gel pad for prophylaxis of pressure sores
ASTOPAD COV150
Heated blanket as under blanket
ASTOPAD COV180
Heated blanket as leg bag



ASTOGEL OPG130
Gel pad for prophylaxis of pressure sores
ASTOPAD COV150
Heated blanket as under blanket
ASTOPAD COV180
Heated blanket as leg bag
ASTOCOVER PAS155
Passive warming insulation
as arm-chest cover

STIHLER ELECTRONIC

Stihler Electronic GmbH

Julius-Hölder-Straße 36 • D-70597 Stuttgart • Germany
www.stihlerelectronic.de • E-Mail: info@stihlerelectronic.de
national: Telefon 07 11- 72 06 70 • Telefax 07 11- 720 67 57
international: Telefon ++ 49- 711- 72 06 70 • Telefax ++ 49- 711- 720 67 57